

DOUBLE DOOR HORIZONTAL AUTOCLAVE

Round & rectangular Double Door Horizontal Steam Sterilizer



Double door horizontal autoclave designed for cleanroom, Operation Theater of Hospital & Microbiology laboratory. This type of autoclave consists of one sterile door and second doors are a non-sterile door. The front door consists of complete control panel inclusive of PID controller, Multiport valve & Pressure gauges where a rear door has a Pressure gauge only. Rear door has no controlling except door opening after pressure comes down to zero.

We design for your Requirement

- Cleanroom
- Microbiology Lab
- Hospitals
- Pharmaceuticals
- Culture & Media sterilization
- Food Industry
- Injectables Manufacturing

Features

- SS 304 construction
- Fully automatic
- Inbuilt Steam generation
- PID controller
- Automatic Purging
- SS 316 Jacket
- Triple Walled Construction
- Radial Locking
- Door interlocking
- Digital Pressure Indicator



Double Door Round Horizontal Autoclave Dimensions

Model No.	Dia. x Height (mm)	Volume (Liters)
Clavex 140D	400 x 1100	140
Clavex 175D	500 x 900	175
Clavex 230D	500 x 1200	230
Clavex 339D	600 x 1200	340
Clavex 552D	750 x 1200	552

Double Door Rectangular Horizontal Autoclave Dimension

Model No.	W x H x D (mm)	Volume (Liters)
Hiclave 182D	450 x 450 x 900	182
Hiclave 243D	450 x 600 x 900	243
Hiclave 324D	600 x 600 x 900	324
Hiclave 432D	600 x 600 x 1200	432
Hiclave 540D	600 x 750 x 1200	540
Hiclave 810D	600 x 900 x 1500	810
Hiclave 1215D	900 x 900 x 1500	1215

