

VERTEX SCIENTIFIC & LAB INSTRUMENTS CO.

(Regd.) : 38/A, Street No.3 Ambedkar Nagar, Haider Pur, Delhi-110088(India)

(Corp.): A223, Ground Floor, Near by Bank of Baroda ATM, Karam Pura Delhi-110015

(Tel.) +91-9999 5737 81 / +91-9999 5737 85 (Email) info@vertexinstruments.com

(Web) www.vertexinstruments.com / www.vertexscientific.com / www.soiltestinginstrument.com

ROCKWELL CUM BRINELL HARDNESS TESTING MACHINE



MAKE: VERTEX | MODEL: VSLIC-RBH.60 | MAKE IN INDIA

Technical Specifications

1. Test Loads (Kgf)	:	60,100,150 (Rockwell) 187.5,250 (Brinell)
2. Initial Loads (kgf)	:	10
3. Max. Test Height(mm)	:	295
4. Depth of Throat(mm)	:	148
5. Max. Depth of Elevating Screw Below Base(mm)	:	295
6. Size of Base (mm) (Approx.)	:	475x220
7. Machine Height (mm)	:	865
8. Net. Weight kg. (Approx.)	:	125

M/s Vertex Scientific & Lab Instruments Co.

Authorised Signatory

VERTEX SCIENTIFIC & LAB INSTRUMENTS CO.

(Regd.) : 38/A, Street No.3 Ambedkar Nagar, Haider Pur, Delhi-110088(India)

(Corp.): A223, Ground Floor, Near by Bank of Baroda ATM, Karam Pura Delhi-110015

(Tel.) +91-9999 5737 81 / +91-9999 5737 85 (Email) info@vertexinstruments.com

(Web) www.vertexinstruments.com / www.vertexscientific.com / www.soiltestinginstrument.com

Standard Accessories

1. Rockwell Diamond Indentor	:	1 No.
2. Superficial Diamond Indentor	:	1 No.
3. Testing Table 50mm	:	1 No.
4. Testing Table 40mm with 'v' Groove for round jobs 6 to 45 mm	:	1 No.
5. Steel Ball Indentor 1/16"	:	1 No.
6. Steel Ball Indentor 2.5 mm	:	1 No.
7. Steel Ball Indentor 5 mm	:	1 No.
8. Test Block HRC	:	1 No.
9. Test Block HRB	:	1 No.
10. Test Block HB 2.5 mm/187.5 Kgf	:	1 No.
11. Test Block Rockwell Superficial 30N	:	1 No.
12. Steel Ball 1/16"	:	5 No's
13. Allen Spanners	:	3 No's
14. Clamping Device	:	1 No.
15. Brinell Microscope	:	1 No.
16. PVC Cover for Machine	:	1 No.
17. Rubber Bellow	:	1 No.
18. Instruction Manual	:	1 Book

(Supply, Installation & Training)

M/s Vertex Scientific & Lab Instruments Co.

Authorised Signatory